




---

# Getter Radon Emanation Measurement

John Harrison  
LZ Lab Technician  
November 2018



# Sample



Valve S2

Getter Cabinet with Portable Panel

Getter cabinet does not heat up unless internal valves are open, so during emanation getter interior was open to the transfer lines, marked in red (up to the valves circled in red). Parts of the transfer lines have been emanated and have emanation rates negligible compared to the Getter. At the end of the second transfer a small seat leak was measured at valve S2 (see slide 11 for details).

Emanated at about 420°C, and with a pressure of 12 psig (roughly 1400 Torr).

Total volume of the transfer lines is 0.6 liters (calculated by inner diameter and length), total internal volume of the getter is 6.8 liters (calculated by volume sharing and associated change in pressure). For volume calculation details see slide 12.




# Overview

---

- Emanated the getter 2 times between August and October, 2018.
- Cabinet and portable panel used for transfer leak-checked using helium leak checker. There were no detectable (external) leaks anywhere in the setup.
- During second transfer pressure dropped by 4.5%. Some portion of this could be lost sample through a valve seat leak (see slide 11 for details).
- Measurements indicated an emanation rate of  $2.26^{+0.28}_{-0.27}$  mBq for all the getter and the cabinet.
  - This is consistent with the preassay assumptions of  $4.69^{+0.35}_{-3.79}$  mBq which was based on Xenon1T radon emanation measurements of two different getters.
    - LZ is estimated to use 3.5x as much getter as Xenon1T.
    - First getter had a rate of 1.34 mBq, resulting in a prediction of 4.7 mBq for LZ.
    - Second getter had a rate of 0.26 mBq, resulting in a prediction of 0.9 mBq for LZ.
    - Higher rate was assumed for central value to be conservative.
- Post-assay background of second transfer was not useable (it was done with a different filter box which turned out to not be functioning correctly).




# First Emanation (August 29 – September 17)






# Second Emanation (September 17–October 16)





# Combined Likelihoods





# Conclusion

---

- Getter emanates at  $2.26^{+0.28}_{-0.27}$  mBq.
  - This includes the getter itself and all components exposed to the flow path in the getter cabinet at operating temperature and 12 psig, as well as the transfer lines connecting it to the portable emanation panel.
  - This is consistent with the preassay assumptions,  $4.69^{+0.35}_{-3.79}$  mBq.





# Backup Slides

---






# Gain Correction Line for First Emanation



- Line used for gain correction shown in green, after gain correction shown on the right.
- Sample run could not be gain corrected based off of Po-210 since the Po-218 signal was large enough to interfere with gain correction.



# Second Emanation with post background run (March 27 – April 9)





# Second Emanation Potentially Lost Sample

---

- During the second emanation, pressure dropped by **4.5%** (from 12 psig to 10.8 psig).
  - This could be due to a combination of a leak and temperature changes in the heated getter and lab air, and/or some interaction with the getter components (adsorption may effect how pressure responds to temperature changes in the getter for example).
  - A seat leak through one of the valves was measured with the helium leak checker at the end of the second emanation, with a rate of about  $1 \times 10^{-8}$  Torr\*liters/sec. If this reading is accurate and constant throughout the emanation the pressure drop due to it should've only been 0.01%.
- If the pressure change is due to temperature rather than a leak, either both the getter cabinet pressure and the room temperature would have to change, or a temperature change more than 4.5% would be needed in one of them. For example:
  - If both changed by the same percent, then getter temperature dropped by roughly 31°C, and room temperature by 13°C.
  - If Getter temperature was constant, then room temperature dropped by roughly 62°C.
  - If room temperature was constant getter temperature would have to have dropped by roughly 37°C
- Since temperature changes this extreme were not noticed, temperature alone does not seem like a good explanation.
- Gas lost early in the emanation would have a smaller effect on emanation efficiency since it would carry less radon out.



# Volume Calculations

---

- Volume of transfer lines and internal getter bypass plumbing calculated based on measured length and internal radius.
- Measurements:
  - Inlet connection length is 105 in, OD is ½ in, ID is 0.375 in.
  - Outlet connection length is 82 in, OD is ½ in, ID is 0.375 in.
  - Getter bypass plumbing is 57 in, OD is ¾ in, ID is 0.625 in.
- Volume in the transfer lines and bypass plumbing (which is open to gas when the valves are opened):

$$\left( 105in * pi \left( \frac{0.375in}{2} \right)^2 + 82in * pi \left( \frac{0.375in}{2} \right)^2 + 57in * pi \left( \frac{0.625in}{2} \right)^2 \right) * \frac{1 \text{ liter}}{61.02 \text{ in}^3} = 0.625 \text{ liters}$$


- Getter internal volume calculated by filling the getter vessel pressure to 15.2 psig and then pumping on the transfer lines until they are approximately at vacuum, and then volume sharing between the getter and transfer lines to get a new pressure of 12.7 psig:

$$\left( 15.2psig + 14.7 \frac{psi}{atm} \right) * V_{Getter} = \left( 12.7psig + 14.7 \frac{psi}{atm} \right) * (V_{Getter} + 0.625 \text{ liters})$$

$$V_{Getter} = \frac{\left( 12.7psig + 14.7 \frac{psi}{atm} \right) * 0.625 \text{ liters}}{15.2psi - 12.7psi} = 6.85 \text{ liters}$$




# First Transfer






# First Transfer






# Cumulative Counts First Transfer





# Resolution First Transfer Backgrounds




ROI efficiency is based on this fit rather than a fit to the sample run since there are more counts and no interference from the Po-218 and Po-214 signals.






# Resolution First Transfer Sample



Resolution shown here is from the Po-214 signal since the Po-218 overlap makes it harder to get resolution from the Po-210 signal.




# Negative Log Likelihoods for First Transfer, Po-214






# Negative Log Likelihoods for First Transfer, Po-218






# Second Transfer with post background






# Second Transfer with post background






# Cumulative Counts, Second Transfer with post background






# Resolution Second Transfer Background





# Resolution Second Transfer Sample




ROI efficiency is based on this fit instead of the background even though there are fewer counts, since the resolution is significantly worse in the background run.






# Negative Log Likelihoods for Second Transfer no post background, Po-214






# Negative Log Likelihoods for Second Transfer no post background, Po-218





# Negative Log Likelihoods for Second Transfer with post background, Po-214





# Negative Log Likelihoods for Second Transfer with post background, Po-218

